

DESCRIPTION

Empty baculovirus particles (BVPs) (empty viral capsids without genomes) can be used to measure polyspecificity of therapeutic antibody candidates in *in vitro* assays. This assesses the viability of an antibody as a therapy early on. They are provided as viral particles in phosphate buffered saline (PBS) at a concentration of $1 \times 10^{13-14}$ pfu/mL. 100 μ L is enough to assay 4 x 96-well plates at a final dilution of 1:100.

BACKGROUND

Some antibodies are able to recognize more than one antigen, which is referred to as polyspecificity (Van Regenmortel, M.H.V., 2014). Therapeutic antibodies have different half-lives in the body that can affect their therapeutic efficiency. One reason for fast clearance from the body is the polyspecificity of the antibody. An ELISA assay using non-specific binding to baculovirus particles (BVPs) of therapeutic antibodies showed correlation with the fast *in vivo* clearance of antibodies (Hötzel *et al.*, 2012). Determining the polyspecificity of an antibody early on during therapeutic antibody development helps with decision making whether to pursue the antibody therapy further.

References:

Van Regenmortel, M.H., "Specificity, polyspecificity, and heterospecificity of antibody-antigen recognition", *J Mol Recognit*, 2014, 27(11): p. 627-39.
Hötzel *et al.*, "A strategy for risk mitigation of antibodies with fast clearance", *mAbs*, 2012, 4:6, 753-760.
Xu *et al.*, "Addressing polyspecificity of antibodies selected from an *in vitro* yeast presentation system: a FACS-based, high-throughput selection and analytical tool.", 2013, *Protein Engineering, Design & Selection*, 26:10, pp. 663-670.

FORMULATION

PBS

PROTOCOL

1. Mix 25 μ L BVP stock with 25 μ L of 50 mM sodium carbonate (pH 9.6) in each well.
2. Incubate 16-24 h at 4°C.
3. Aspirate the unbound BVPs from the wells and perform the rest of the steps at room temperature.
4. Add 100 μ L of blocking buffer (PBS with 0.5% BSA) and incubate for 1 hr.
5. Wash 3 times with 100 μ L of PBS.

6. Add 50 μ L of 1 μ M primary antibody and incubate for 1 hr.
7. Wash 6 times with 100 μ L of PBS.
8. Add 50 μ L of anti-Human-IgG-HRP and incubate for 1 hr.
9. Wash 6 times with 100 μ L of PBS.
10. Add 50 μ L of TMB substrate and incubate for 10-15 min.
11. Stop the reaction with 50 μ L of 2M sulfuric acid.
12. Read absorbance at 450 nm.
13. BVP score is determined by normalizing absorbance by control wells with no test antibody.

STORAGE

-80°C